NAME: CHEMISTRY PAPER 3 FORM THR TIME: 1 ¾ HOURS			CLASS:
 You are provided with: Dilute hydrochloric acid labeled Solution B containing 3.15g of solution. Sodium hydroxide solution labeled 	dibasic acid H ₂ C ₂ O ₄	4.2H ₂ O dissolved to	form 500cm ³ of
You are required to: - Standardize the sodium hydro - Use the standardized solution (oncentration of soluti	ion A.
Procedure: (i) Place solution B in a clean burette. (ii) Using a pipette and a pipette filler, por phenolphthalein indicator. Titrate with the solution of pipette used =	th solution B and rec		
	I	II	III
Final burette reading (cm ³)			
Initial burette reading (cm ³)			
Volume of solution B used (cm ³)			
(a) Calculate the:	B used.		(4 mks) (1 mk)
(ii) Concentration of the dibasic	solution B in moles	per litre. (H = 1, O =	= 16, C =12) (4 mks)
(ii) Concentration of the sodium	hydroxide solution	C in moles.	(4 mks)

	Using a clean measuring cylinder, measure 25cm ³ of solution A and transfer into 250ml volumetric flask.
(ii)	Use distilled water to top up the solution in the volumetric flask to the mark.
(iii)	Mix the solution well and label it a solution K.
	Clean the burette and place solution K.
(v)	Using a pipette and pipette filler, place 25cm ³ of solution C into a conical flask and add two drops of phenolphthalein indicator.

(vi) Titrate the concentration in the conical flask with solution K. Record the results in table 2. Repeat

the titration two more times and record your results. $\overline{\textbf{TABLE 2}}$

Volume of pipette used = ______(1 mk)

(c) Determine the concentration of solution A in moles per litre.

	I	II	III
Initial burette reading (cm ³)			
Final burette reading			
Volume of solution K used (cm ³)			

	ar burette reading			
Vo	lume of solution K used (cm ³)			
(a) What is the average volume of solution K used?		(1 mk)		
(b)	Calculate the concentration of solu	ution K in moles pe	r litre.	(4 mks)

(3 mks)

 2. You are provided with solid D. Carry out the tests in the spaces provided. (a) Place solid D in a boiling tube. Add about 8cm³ solution into four portions. 	•
Observations	Inferences
(b) To the first portion, add aqeous sodium hydroxid	e dropwise until in excess. (3 mks)
Observations	Inferences
(c) To the second portion, add aqeous ammonia drop	
Observations	Inferences
(d) To the third portion, add barium chloride solution	
Observations	Inferences

(e) To the fourth portion, add acidified lead (II) nitra	ate solution. (2 mks)
Observations	Inferences