

CONFIDENTIAL
TERM 2 2017
BIOLOGY FORM 3 PAPER 3

Each student should have;

1. 4 test tubes and a test tube rack
2. Iodine solution – supplied with dropper
3. 10cm visking tubing labeled J
4. Piece of string 20cm long
5. 10cm³ solution of a mixture of soluble starch and glucose labeled K

N/B. 30g glucose mixed with 3g starch then add 100cm³ water and heat to boil then cool.

6. 500 ml beaker
7. Adequate distilled water / clean rain water
8. Benedicts solution
9. Means of heating/ Bunsen burner
10. Measuring cylinder – 10 ml

END OF PRINTED PAGE
